




Focus on Delivery Systems

ASI Medical
(800) 566-9953
asimedical.net

BOYD Industries Inc
(800) 255-2693
boydindustries.com

Dexta Corp
(800) 733-3982
dexta.com




Product Name	2021 Standard Mobile System	CSU 356 FLM	SpaceSaver
			
Product type	Self-contained delivery system	Delivery system	Cabinet
Product use	Hygiene, orthodontics, and other two-handed dentistry	Chairside delivery system with instrumentation	Delivery system, storage, computer setup, and countertop
Dimensions	H 35.5"; W 19.5"; D 22" (W and D measurements may vary depending on options)	W 22 3/4"; H 31"	W 20" x D 28" x H 31", including casters
Sold individually or as a set?	Individually	Individually	As a set
Mobile or fixed?	Mobile	Mobile	Mobile
Materials available	N/A	Cabinet: Wilsonart, Nevamar. Formica; Tops: Corian, Silestone, Wilsonart, Nevamar, Formica	Laminates and solid surface tops
Cost	\$8,250	Retail base price \$4,295. Doctor Direct pricing available.	\$2,890 list plus options
Length and cost of warranty	Included in price. ASI warrants products against defects in material and/or workmanship for 1 year from the time of delivery for parts and 90 days for labor.	3 year parts and 1 year labor are included.	3 years parts and 1 year service
Special features	Internal air compressor system, suction system with separate vacuum pump, electric canister purge, aseptic dual closed water bottles, one pneumatic handpiece connection, air/water syringe, titanium saliva ejector, and straight high-volume evacuator are included. The WilsonArt™ solid surface worktop provides an easy-to-clean, aseptic surface.	ETL and FDAA approved, two duplex outlets to accommodate equipment, self-closing maple drawers, undermount guides, superior workmanship and quality. All plywood construction, foot control, regulators, two-handpiece automatic delivery system.	Dexta cabinets come standard with intake and exhaust vents, a thermostat-controlled fan, and a surge suppressor. Most models can even accept a full-size CPU.

Focus on Delivery Systems

DOME Inc
(800) 423-5767
domeortho.com

J. Morita USA
(877) JMORITA (566-7482)
morita.com/usa

Summit Dental Systems
(800) 275-3368
summitdental.com

Product Name	Coronado	Spaceline EMCIA	900SH 2HP Mobile Cabinet
			
Product type	Operatory delivery unit	Delivery system	Mobile delivery cabinet
Product use	Provides instrumentation, storage and computer workstation	Chair and delivery system	Delivery unit for orthodontic use
Dimensions	W 18" to 26.5"	Approx L 72" x W 26"	H 32 3/4" x W 25 3/8" x D 18 1/2"
Sold individually or as a set?	Individually	As a set. System may be customized based on individual preference.	Individually
Mobile or fixed?	Mobile	Fixed	Mobile
Materials available	Corian, Wilsonart, Formica, Nevamar, and others	Nine colors to choose from	Custom laminate or solid surface available
Cost	Varies; many options available	Suggested retail \$35,500	\$2,780 MSRP
Length and cost of warranty	5-year limited warranty	Free 1-year warranty	5 years on the delivery unit and 1 year on the cabinet
Special features	Outstanding quality and design	The Spaceline EMCIA was designed to promote a balanced posture, while reducing operator stress and fatigue. The handpieces are positioned on the back of the chair below the headrest, while soft-touch switches are located on both the backrest and tray to help decrease unnecessary movements and body contortion. Other features include warm air/water, a low chair height setting (approximately 18"), cuspidor, and piezo electric scaler.	Cabinets are available in a standard white color with drawers that can be painted in white, gray, surf, or almond. They can also be furnished with any laminate option of your choosing for an additional fee. This laminate can be applied to either the whole cabinet or to individual sides, doors, or drawers of the cabinet. The option to add a solid-surface top is also available for an additional fee.