


Focus On RETAINERS

Accutech Orthodontic Lab Inc
accutechortho.com
(800) 734-7855

Align Technology
aligntech.com
(408) 470-1000

AOA Labs
aoalab.com
(800) 262-5221

BioMers LLC
simpliclear.com
(904) 519-7778




Product name	QuickBond BLR (Bonded Lingual Retainer)	Vivera Retainers	Retainers	ASTICS Retainer
Type of retainer	Bonded	Clear plastic	Active/spring, Bonded, Clear plastic, Hawley	Hawley
Material	Composite, Wire	Plastic	Acrylic, Composite, Plastic, Wire	Acrylic, Composite, Wire
Turnaround time	3 days in lab	Varies	Varies by design	N/A
Initial cost	\$46.40	Varies	Costs vary per design. Call for specifics.	\$12 per unit x 100 units (minimum order quantity) = \$1,200
Repair or remake cost	Zero	Varies	Varies depending on warranty. Call for specifics.	N/A
How long is it designed to last?	Until patient wants it removed	Varies	5 years	2 years
Length of warranty	90 days	Varies	3 months	N/A
Aesthetic options	Aesthetic wire	N/A	X	X
	Clear plastic	N/A	X	
	Labeling	N/A		
	Lingual wire placement	N/A		X
	Multiple colors	N/A		X
	Personalized images	N/A		X
Special features	This appliance comes with custom composite pads on each tooth and a custom dual-tray delivery system.	Can be set up on one-time or annual shipments.	All retainer images can be fabricated from digital impressions/scans provided by your office.	Designed to work similar to traditional metal Hawley retainers. Durable, effective, gentle, and nonabrasive. It can be fabricated to the following designs according to doctors' preferences: Wrap-around and Soldered.

ClearCorrect
clearcorrect.com
(888) 331-3323

DynaFlex
dynaflex.com
(800) 489-4020

**DENTSPLY
Raintree Essix**
essix.com
(800) 883-8733

**Forestadent
USA**
forestadentusa.com
(800) 721-4940



Product name	ClearCorrect Retainer	QCM Retainer	MTM® Clear•Aligner	TRACK A Thermoforming Foils
Type of retainer	Clear plastic	Hawley	Clear plastic	Clear plastic
Material	Plastic	Acrylic, Plastic	Plastic	Plastic
Turnaround time	Less than 1 week	1 week in laboratory	N/A (made in office)	Less than 15 minutes
Initial cost	\$50 per arch	\$79.60	N/A	Please call for pricing
Repair or remake cost	\$50 per arch	Acrylic repair \$20	Less than \$10	Cost of one Thermofoil
How long is it designed to last?	Three months when used as recommended. Can vary based on doctor-prescribed patient use and patient wear times.	2 years full-time wear	18 to 24 months (varies based on selected plastic)	With proper care, designed to last the life of the retainer.
Length of warranty	90 days	6 months	N/A	Limited warranty. Call for details.
Aesthetic options	Aesthetic wire			
	Clear plastic	X	X	X
	Labeling		X	
	Linqual wire placement			
	Multiple colors		X	
	Personalized images		X	
Special features	ClearCorrect retainers are made with Zendura®. Dental models (the production molds) are included at no additional charge so doctors can easily make replacements in-house if the retainers are lost or wear out.	Crystal-clear labial plastic replaces traditional stainless steel labial bar.	MTM® Clear•Aligner In-Office utilizes Essix® Brand plastics and the Hilliard Thermoplier® Pliers to fabricate clear plastic aligners that move teeth with integrated "force points."	Do not require drying or bond-enhancing mediums. Each foil has an integrated separating foil, reducing time and material during preparation. Available in five dimensions. Used for invisible and standard retainers, and splints.

Focus On RETAINERS

**G&H
Orthodontics**
ghortho.com
(800) 526-1026

**Gold'n Braces
Inc**
goldnbraces.com
(727) 785-5491

**Great Lakes
Orthodontics Ltd**
greatlakesortho.com
(716) 871-1161

**Les Idées
LLC**
lesideesinc.com
(760) 650-6570



Product name	everstick® Ortho	ZMRW	Hawley Retainer	Retain-It™ (retainer storage)
Type of retainer	Bonded	Bonded	Hawley	N/A
Material	Composite	Wire	Acrylic, Wire	N/A
Turnaround time	Chairside application	N/A	5 days	7 to 10 days
Initial cost	Contact company	N/A	\$62.05	\$100 for 10-piece pack
Repair or remake cost	N/A	N/A	Acrylic: No charge. Wire: Under 90 days. No charge.	N/A
How long is it designed to last?	Permanent	N/A	Many years if properly taken care of	2 years with specific components longer
Length of warranty	N/A	N/A	Acrylic: Lifetime. Wire: 90 days.	30 days
Aesthetic options	Aesthetic wire	N/A	X	X
	Clear plastic	N/A		X
	Labeling	N/A		X
	Lingual wire placement	N/A	X	
	Multiple colors	N/A		X
	Personalized images	N/A		X
Special features	everstick® Ortho incorporates a thousand glass fiber strands bound together with a preimpregnated light cure adhesive that allows adaptation directly in the patient's mouth. The integrated fiber bundle is both malleable and tacky. .030" (.75 mm) passive diameter.	24-carat gold-plated retainer wire	Many different patterns and colors to choose from.	Made with machine-washable neoprene, colored zippers, an elastic band to hold the retainer case in place, and a heavy-duty carabineer for universal attachment. This product reduces loss and damage of retainers, simplifies storage, and increases patient compliance.

Focus On RETAINERS

Ortho Organizers Inc

orthoorganizers.com
(800) 547-2000

PAR Orthodontics Lab

parortholab.com
(949) 472-4788

Pilot Orthodontics Inc

pilotortho.com
(800) 352-0949

Ribbon Inc

ribbon.com
(800) 624-4554



Product name	Hawley Spring Retainer	Schwartz	Modified Spring Aligner with two-step design	Ribbon Reinforcement Ribbon	
Type of retainer	Active/spring	Hawley	Active/spring	Bonded	
Material	Plastic, Wire	Acrylic, Wire	Acrylic, Wire	Composite	
Turnaround time	5 days	3 days in lab	3 days	Made direct	
Initial cost	N/A	\$64.50	\$61.50 plus \$5.50 each reset tooth	\$290.69 for 20 3-3 retainers	
Repair or remake cost	N/A	No charge guarantee for 90 days if faulty or defective	Our fault: free. Otherwise, \$30 minimum.	Contact company	
How long is it designed to last?	Depends on patient care	With good care will last a year or more	1 year past alignment goal	Longevity is generally determined by bond of composite to teeth	
Length of warranty	N/A	90 days	N/A	Contact company	
Aesthetic options	Aesthetic wire			N/A	
	Clear plastic	X		X	N/A
	Labeling	X		X	N/A
	Lingual wire placement				N/A
	Multiple colors	X	X	X	N/A
	Personalized images	X	X	X	N/A
Special features	Free retainer case with every retainer.	Ideal for correcting cross-bites, and helps eliminate anterior crowding. Available for upper and lower arch. Posterior bite blocks are an option.	Patient's name in all removable appliances. Choice of clasping. No charge for solid colors (38 choices) and glitter (11 choices). 1-, 2-, 3-, and 4-step designs available for additional cost.	Directly made fiber-reinforced composite resin retainers.	

Focus On RETAINERS

**SML Space
Maintainers Lab**
SMLglobal.com
(800) 423-3270

**Specialty Appliances
Orthodontic Laboratory**
specialtyappliances.com
(800) LAB-INFO (522-4636)

**Tortorici Orthodontic
Lab Inc**
ClearFusionOrtho.com
(631) 642-1900

**TP
Orthodontics**
TPOrtho.com
(800) 348-8856



Product name	Inman Aligner	Fixed Lingual Retainer	Clear Fusion Hybrid Retainer	Spring Aligner Plus
Type of retainer	Active/spring	Bonded	Clear plastic	Active/spring
Material	Acrylic, Wire	Composite, Wire	Acrylic, Plastic	Acrylic, Wire
Turnaround time	5 days	3 days in the lab	3 days in-house	1 week
Initial cost	\$189	\$38 to \$52	\$40 to \$65 per arch	\$85
Repair or remake cost	Contact company	No charge for 1 year	N/A	\$28.50 for repair or \$84.80 for remake
How long is it designed to last?	Indefinitely	Indefinitely	1 year plus with proper care	Indefinitely
Length of warranty	91 days	1 year	90 days	For length of treatment. Does not cover patient neglect/abuse.
Aesthetic options	Aesthetic wire	X	X	
	Clear plastic			X
	Labeling	X		X
	Lingual wire placement		X	X
	Multiple colors	X		X
	Personalized images	X		
Special features	Nickel titanium coil springs power two aligner bows which gently oppose each other, guiding the teeth into their new position. Most cases are completed within 6 to 18 weeks.	The appliance can be constructed with scalloped wire to prevent flossing issues.	New generation of wireless, clear, and aesthetic appliances with rigid lingual support. Clear Fusion Appliances can be made in many different rainbow colors or clear.	The Palatal/Lingual portion is segmented, allowing the appliance to seat positively because of the flexibility of the lingual pad on the incisors without any compromise to appliance strength. Customers can achieve flawless alignment and correction, especially involving labial movement.

Tru-Tain Orthodontic & Dental Supply
 tru-tain.com
 (800) 843-0904



Tru-Tain DX Retainer Material

Clear plastic

Plastic

Immediately made in office

\$1.15

\$1.15

6 months to 1 year


N/A

X

It can be made in-office immediately after debonding orthodontic brackets.

Receive your
FREE
 subscription to

Orthodontic Products connects orthodontists and dentists who perform orthodontic procedures to product manufacturers.



22 Finishing Archwires
 30 Interviewing Techniques
 36 Focus on Functional Appliances

SCAN THE QR CODE TO SUBSCRIBE OR RENEW!
 (<http://antheimedia.mkt5097.com/OPReg/Registration/>)

ORTHODONTIC PRODUCTS

SCAN ME!



Scan the QR code above with your smartphone to connect to our Market Guide.

ORTHODONTIC PRODUCTS

“The most efficient and predictable Class III corrector appliance. It will change the way you diagnose and treat Class III patients.”

Dr. Dana Fender, DMD, MS
 FENDER-GOGGANS ORTHODONTICS

CS-2000®
 Your Class III Solution



Patented Archwire Attachment

CS-2000® 5 patient kit

Easy To Place

Little To No Maintenance

Up to 5mm of Correction

350 Grams of Low Continuous Force

Proven Results!



Case photos compliments of Dr. Andy Hayes of St. Louis, MO

DynaFlex

www.dynaflex.com | 800-489-4020

060413 © 2013 DynaFlex®, St. Louis, MO 63074. All rights reserved.
 CS-2000® patent #6719557 Pivot with slotted base patent #6036488